

Sleeper dB Blocker®

The Sleeper dB Blocker[®] can solve one of the important health issues facing many people today – lack of deep, uninterrupted sleep.

$$\label{eq:second} \begin{split} & \mathbf{I} \text{ is custom fit to the individual wearer for maximum comfort using Skinsoft}^{\mathsf{TM}} \\ & \text{medical grade silicone. The Sleeper dB Blocker}^{\circledast}, when inserted allows wearer to place head on pillow without causing pressure points. Because of custom molding and design, it will not come out during sleep. The wearers can still hear alarm clocks, fire alarms, smoke detectors and buzzers. Also available in a vented model for situations where you need to be discreet. \end{split}$$

RECOMMENDED USES

- Those with problems sleeping
- Shift Workers
- Persons with snoring spouses
- Persons living near high traffic areas
- Flight Attendants
- Bartenders/Servers

The use of this hearing protector is not limited to the above applications. Where noise protection and discretion are required the Sleeper dB Blocker[®] is recommended.

FEATURES

- Smallest hearing protector available
- Blocks out snoring
- Inconspicuous in the ear
- Available in both solid and vented which enables communication



Comfortable to lie on

<u>BENEFITS</u>



- Ideal for shift workers because it allows wearer to have an undisturbed nights sleep which makes them safer on the job
- Spouses are able to sleep through the night without annoyance
- Persons who work with the general public such as servers, bartenders, flight attendants, and emergency response personnel can protect their hearing without appearing to shut out the public

<u>SPECIFICATIONS</u>

- Made of SkinsoftTM medical grade silicone
- Hypoallergenic
- Non-porous
- Durometer of 30 (800 PSI tear strength)
- 25 lbs tear strength on metal anchor
- Product inert to 500 F
- Smallest HPD available for comfort when lying down or discretion in public jobs where hearing protection is required
- Colours Available Blue, Green, Orange, Red, Black, Purple, and Clear

CARE AND MAINTENANCE

- Lubrication with Petroleum Jelly essential for first 2 weeks
- Wash regularly with soap and water
- Re-lubricate as needed after washing

